

THIN TIPS FOR HIGHEST PRECISION

FREE OF
detectable
DNase

FREE OF
detectable
human DNA

FREE OF
detectable
RNase

 non-
pyrogenic

 PCR

GEL-LOAD PIPETTE TIPS


Pipetting liquids into narrow vessels or for loading gels requires high precision and a steady hand. To facilitate this process, Greiner Bio-One offers gel-load pipette tips with a slim shape that allows even the thinnest gels to be loaded without damaging them. Due to the long and flexible tip, pipetting at the bottom of very narrow vessels is also possible without any problems.

The gel-load pipette tips are packed in a resealable bag, which means that the tips are safely protected from possible contamination even after the bag has been opened.

- / Volume range: 1 - 200 μ l
- / Resealable bags with 1,000 tips each
- / Free of detectable DNase, RNase, human DNA, pyrogens and PCR inhibitors

KEY FACTS

- / Made of polypropylene
- / Universal fit for all common pipettes
- / 5 μ l graduation for easy control of the pipetted volume
- / Tip with 0.57 mm outer diameter suitable for loading gel pockets
- / 83 mm total length to reach even the bottom of narrow vessels



PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERING
PLEASE VISIT OUR WEBSITE OR CONTACT US.

Gel-Load Pipette Tips

Volume: 200 μ l, Volume range: 1 - 200 μ l, Type of packaging: bag

Item no.	Product colour	Qty. inner / outer
775390	<input type="radio"/> natural	1.000 / 10.000