



TRIPLE SAFETY FOR YOUR CLEANROOM

STERILE	FREE OF detectable DNase	FREE OF detectable human DNA	FREE OF detectable RNase
	non-cytotoxic	non-pyrogenic	


TRIPLE PACKED PRODUCTS

Triple packed products are used where sterile transfer of products is essential. With single packaging, contamination adhering to the outside of the bag can be carried into clean room areas. Therefore all single packaging materials must be wiped clean with 70 % alcohol or other disinfectants. Triple packed products eliminate the need for such cleaning procedures and prevent contamination from reaching high sensitive production areas.

- / Convenient and easy clean room handling
- / Facilitates workflows
- / Reduced processing time
- / Three individual peel bags
- / High sterilisation level SAL 10⁻⁶

PACKED IN THREE INDIVIDUAL PEEL BAGS

- / Easy-to-open peel bags
- / No need for outer cleaning procedures
- / Increasing bag size for safer handling
- / Inner bag is labeled with product traceability information



PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERING
PLEASE VISIT OUR WEBSITE OR CONTACT US.

Triple-packed

CELLSTAR® Standard Cell Culture Flasks

Raw material: PS, Surface treatment: TC, Triple-packed: yes, Sterile: +

Item no.	Flask design	Growth area	Cap colour	Working volume	Total volume	Cap design	Sterile	Qty. inner / outer
690160-TRI		25 cm ²	● red	5 ml - 10 ml	50 ml	screw cap	+	10 / 130
690175-TRI		25 cm ²	● red	5 ml - 10 ml	50 ml	filter screw cap	+	10 / 130
658170-TRI		75 cm ²	● red	15 ml - 38 ml	250 ml	screw cap	+	5 / 90
658175-TRI		75 cm ²	● red	15 ml - 38 ml	250 ml	filter screw cap	+	5 / 90
660160-TRI	flat	175 cm ²	● red	20 ml - 45 ml	550 ml	screw cap	+	5 / 40
660175-TRI	flat	175 cm ²	● red	20 ml - 45 ml	550 ml	filter screw cap	+	5 / 40

Triple-packed

Suspension Culture Flasks

Raw material: PS, Surface treatment: suspension, Triple-packed: yes, Sterile: +

Item no.	Flask design	Cap colour	Total volume	Cap design	Sterile	Qty. inner / outer
690190-TRI		○ white	50 ml	screw cap	+	10 / 130
690195-TRI		○ white	50 ml	filter screw cap	+	10 / 130
658190-TRI		○ white	250 ml	screw cap	+	5 / 90
658195-TRI		○ white	250 ml	filter screw cap	+	5 / 90
660190-TRI	flat	○ white	550 ml	screw cap	+	5 / 40
661195-TRI	high	○ white	650 ml	filter screw cap	+	4 / 28

Triple-packed

Cell Culture Dishes

Surface treatment: TC, Vent nock: yes, Triple-packed: yes, Sterile: +

Item no.	Height	Ø nominal size	Growth area	Working volume	Total volume	Sterile	Qty. inner / outer
628160-TRI	15 mm	60 mm	21 cm ²	6 ml - 7 ml	28 ml	+	10 / 300
664160-TRI	20 mm	100 mm	58 cm ²	16 ml - 17 ml	100 ml	+	15 / 180
639160-TRI	20 mm	145 mm	143 cm ²	25 ml - 27 ml	240 ml	+	5 / 120

Triple-packed

CELLSTAR® Cell Culture Tubes

Support skirt: no, Graduation: yes, Writing field: yes, Raw material: PP, Bottom shape: conical, Cap design: screw cap, Type of packaging: triple packed, Sterile: +

Item no.	Height	Ø	Product colour	Cap colour	Colour writing field	Graduation colour	Working volume	Nominal volume	Sterile	Qty. inner / outer
188271-TRI	120 mm	17 mm	○ natural	● blue	○ white	● blue	≤14 ml	15 ml	+	5 / 500
227261-TRI	115 mm	30 mm	○ natural	● blue	○ white	● blue	≤50 ml	50 ml	+	5 / 250

PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERING
PLEASE VISIT OUR WEBSITE OR CONTACT US.

Triple-packed Serological Pipettes

Raw material: PS, Packaging: plastic / plastic, Triple-packed: yes, Sterile: +

Item no.	Description	Feature	Graduation	Nominal volume	Packaging	Triple-packed	Sterile	Qty. inner / outer
710183-TRI	2 ml aspiration pipette	no plug	-	2 ml	paper / plastic	yes	+	10 / 100
604160-TRI	1 ml pipette		1/100	1 ml	plastic / plastic	yes	+	10 / 100
710160-TRI	2 ml pipette		1/100	2 ml	plastic / plastic	yes	+	10 / 100
606160-TRI	5 ml pipette		1/10	5 ml	plastic / plastic	yes	+	10 / 100
607160-TRI	10 ml pipette		1/10	10 ml	plastic / plastic	yes	+	10 / 100
760160-TRI	25 ml pipette		2/10	25 ml	plastic / plastic	yes	+	10 / 100
768160-TRI	50 ml pipette		1/2	50 ml	plastic / plastic	yes	+	10 / 100

Triple-packed Cryos

Starfoot: yes, Graduation: yes, Writing field: yes, Raw material: PP, Cap design: screw cap, Triple-packed: yes, Sterile: +

Item no.	Height	Ø	Barcode	Bottom shape	Cap colour	Working volume	Thread type	Sterile	Qty. inner / outer
123263-TRI	42 mm	12.5 mm	no	conical	○ natural	≤1.2 ml	internal	+	10 / 200
122263-TRI	48 mm	12.5 mm	no	round	○ natural	≤2 ml	internal	+	10 / 200
126263-TRI	47 mm	12.4 mm	no	round	○ natural	≤2.2 ml	external	+	10 / 200