

SPECIAL
DISCOUNT*


greiner
BIO-ONE

SECURE THE FUTURE OF YOUR CRYO SAMPLE STORE

Triple Coded Cryo.s Tubes

Effortless sample identification and seamless upscaling

Explore our triple-coded Cryo.s cryogenic storage tubes! Tailored for your needs, these tubes support automation, ensure sterility, and offer flexible working volumes. Trust in our CE IVDR-certified and airfreight-ready solutions to advance your cryopreservation processes.

- / Eliminate errors due to illegible handwriting and tube labeling
- / Prevent sample mix-ups with unique IDs
- / Start tracking instantly with human-readable tube IDs
- / Effortlessly scale your throughput with machine-readable codes

KEY FACTS:

- / Available in 1, 2, and 4 ml working volumes
- / Unique, serialized tube ID
- / Human-readable ID format
- / Standard or customized barcodes
- / Additional ID representation through linear barcodes and Datamatrix ECC200 codes
- / Suited for LN2 vapour phase storage (-196°)
- / Made from USP Class VI, medical-grade materials
- / Sterility assurance 10⁻⁶
- / CE IVDR certified
- / Airfreight certified
- / Automation compatible



REQUEST YOUR INDIVIDUAL OFFER
OR SAMPLE FOR TESTING



Greiner Bio-One GmbH Frickenhausen, Germany
PHONE +49 7022 948-0 / FAX +49 7022 948-514 / E-MAIL info@de.gbo.com
Greiner Bio-One is a global player.
Find the contact details of your local partner on our website.

www.gbo.com

PRODUCT OVERVIEW

REQUEST YOUR INDIVIDUAL OFFER
OR SAMPLE FOR TESTING

Triple Coded Cryo.s Tubes

Starfoot: yes, Barcode type: linear barcode and datamatrix code,
Raw material: PP, Cap design: screw cap, Sterile: +, Cap colour: natural

Item no.	Height	Ø	Bottom shape	Working volume	Thread type	Qty. inner / outer
123263-2DG	42 mm	12.5 mm	conical	≤1.2 ml	internal	100 / 500
122263-2DG	48 mm	12.5 mm	round	≤2 ml	internal	100 / 500
126263-2DG	47 mm	12.4 mm	round	≤2.2 ml	external	100 / 500
127263-2DG	83 mm	12.4 mm	round	≤4 ml	external	50 / 300



GET YOUR INDIVIDUAL DISCOUNT AND/OR SAMPLE TODAY!*

* This discount is available for new customers on their first order of this product in specific regions.

For inquiries, please use the contact form accessible via the QR code or link.