



**FOR
MAN'S
BEST FRIEND**



MINICOLLECT®

Veterinary tubes

MiniCollect® tubes for veterinary use can be used with large, small and medium-sized animals. The MiniCollect® tube has an integrated scoop and enables drops of blood to be collected easily and hygienically.

The tube's cap can be pierced by an analysis needle while still closed and remaining leakproof. It is not absolutely necessary to remove the cap first. This ensures hygienic processing of the samples.

- / In combination with carrier tubes (13 x 75 mm), the tubes can be processed in a standard centrifuge.
- / Label text is colour coded so that the tube type can be identified even when the cap has been removed.
- / Transparent label for optimal visual inspection of the sample material.

MINICOLLECT® COMPLETE WITH PRE-ASSEMBLED CARRIER TUBES.

MiniCollect® Complete tubes are CLSI AUTO1-A-compliant. The tubes fit into 13 x 75 mm racks and are suitable for automatic processing in the laboratory.



PRODUCT OVERVIEW

FOR FURTHER INFORMATION OR TO ORDER SAMPLES
PLEASE VISIT OUR WEBSITE OR CONTACT US.

MiniCollect® Complete

for veterinary use only

Item No.	Volume range	Cap colour	Additive	Tube size	Label type	Inner / Outer [Qty.]
450547VET	0.25 / 0.5ml	● lavender	K2E K2EDTA	13 x 75	transparent	500 / 2,000
450555VET	/	○ white	Z No Additive	13 x 75	transparent	500 / 2,000

MiniCollect® Tubes

for veterinary use only

Item No.	Volume range	Cap colour	Additive	Tube size	Label type	Inner / Outer [Qty.]
450533VET	0.5 / 0.8ml	● gold	CAT Serum Sep Clot Activator	10 x 43	transparent	500 / 2,000
450571VET	0.5 / 0.8ml	● gold	CAT Serum Sep Clot Activator	10 x 43	paper	500 / 2,000
450535VET	0.8ml	● light-green	LH Lithium Heparin Separator	10 x 43	transparent	500 / 2,000
450572VET	0.8ml	● light-green	LH Lithium Heparin Separator	10 x 43	paper	500 / 2,000
450537VET	1ml	● green	LH Lithium Heparin	10 x 43	transparent	500 / 2,000
450574VET	1ml	● green	LH Lithium Heparin	10 x 43	paper	500 / 2,000
450532VET	0.25 / 0.5ml	● lavender	K2E K2EDTA	10 x 43	transparent	500 / 2,000
450573VET	0.25 / 0.5ml	● lavender	K2E K2EDTA	10 x 43	paper	500 / 2,000
450542VET	/	○ white	Z No Additive	10 x 43	transparent	500 / 2,000

MiniCollect® Carrier Tubes

Item No.	Colour code	Thread type	Tube size	Inner / Outer [Qty.]
450417	○ transparent	PREMIUM	13 x 75	100 / 500
450418	● amber	PREMIUM	13 x 75	100 / 500

MiniCollect® Extra Cap

Item No.	Colour code	Inner / Outer [Qty.]
370267B	○ white	1 / 9,008