

# ENHANCED SAFETY FOR PATIENTS AND CLINICIANS

## VACUETTE® SAFETY Blood Collection Sets

This product provides a high degree of safety as the safety mechanism is activated while the cannula is still inside the vein. Proper activation by the user is associated with a clearly audible click. The SAFETY Blood Collection Set thereby provides enhanced protection against needlestick injuries.

The SAFETY Blood Collection Set is a sterile, single-use winged needle, connected to flexible tubing, with or without a luer adapter, or preassembled with a standard or blood culture holder.

- / A view area in the body of the device provides a visual indication of successful venipuncture.
- / Tubing availability in different lengths offers more flexibility in blood collection.
- / The safety mechanism provides simple in-vein activation.
- / May be used for venous blood collection or infusions for up to five hours.\*

### ERGONOMIC PACKAGING IS EASY TO OPEN AND MINIMIZES WASTE.

Sustainability matters to Greiner Bio-One. Thoughtfully designed packaging for our SAFETY Blood Collection Sets optimizes the units per shelf box, case and pallet, which not only reduces our carbon footprint but also offers logistical advantages for our customers.

### THE NUMBERS SPEAK FOR THEMSELVES:

- / 45.39% fewer single-use pallets required for shipping per year
- / 74 fewer container trucks per year
- / 286.58 lbs reduction in packaging material (shelf box and shipping carton) per year



## PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERS  
PLEASE CONTACT US OR VISIT OUR WEBSITE  
[www.gbo.com](http://www.gbo.com)

### SAFETY Blood Collection/Infusion Set

Item No.	Color Code	Needle Size	Needle Length	Needle Wall Type	Tubing Length	Bag [pcs]	Shipping carton [pcs]
450091V1	● green	21 G	0.75"	thin wall	7.5"	90	1,080
450092V1	● blue	23 G	0.75"	thin wall	7.5"	90	1,080
450090V1	● orange	25 G	0.75"	thin wall	7.5"	90	1,080

### SAFETY Blood Collection Set + Luer Adapter

Item No.	Color Code	Needle Size	Needle Length	Needle Wall Type	Tubing Length	Bag [pcs]	Shipping carton [pcs]
450081V1	● green	21 G	0.75"	thin wall	7.5"	90	1,080
450095V1	● green	21 G	0.75"	thin wall	12"	90	1,080
450084V1	● blue	23 G	0.75"	thin wall	4"	90	1,080
450082V1	● blue	23 G	0.75"	thin wall	7.5"	90	1,080
450096V1	● blue	23 G	0.75"	thin wall	12"	90	1,080
450099V1	● orange	25 G	0.75"	thin wall	12"	90	1,080

### SAFETY Blood Collection Set + Holder

Item No.	Color Code	Needle Size	Needle Length	Needle Wall Type	Tubing Length	Bag [pcs]	Shipping carton [pcs]
450085V1	● green	21 G	0.75"	thin wall	7.5"	35	420
450160V1	● green	21 G	0.75"	thin wall	12"	35	420
450086V1	● blue	23 G	0.75"	thin wall	7.5"	35	420
450161V1	● blue	23 G	0.75"	thin wall	12"	35	420

### SAFETY Blood Collection Set + Blood Culture Holder

Item No.	Color Code	Needle Size	Needle Length	Needle Wall Type	Tubing Length	Bag [pcs]	Shipping carton [pcs]
450184	● green	21 G	0.75"	thin wall	12"	24	240
450185	● blue	23 G	0.75"	thin wall	12"	24	240

\* Except SAFETY Blood Collection Set + Blood Culture Holder.