

Stonehouse, May 2022

INFOMAILING

PRODUCT MODIFICATION: CRYO.S WITH DATAMATRIX CODE AND LINEAR BARCODE / 2 ML VERSIONS WITH INTERNAL THREAD

Dear Sir or Madam,

We are continuously striving to optimise our product portfolio and to improve existing products. In this context we would like to inform you about a design change of the **Cryo.s with Datamatrix Code and Linear Barcode / 2 ml versions with internal thread**.

In the future these products will have the following properties:

- Produced with next generation manufacturing technology utilising two-component injection moulding in combination with laser-labelling
- Optimum stability and robustness of tube-marking due to 100% integration of barcoding areas into tube base body
- Unambiguous tube coding with Datamatrix code at bottom and corresponding linear barcode with human readable feature on side
- Graduation with 0.5 / 1.0 / 1.5 and 2.0 ml-marks
- Unchanged functionality and automation compatibility

[more >](#)



Fig. left: Current Design



**Fig. right: Future design
(details may be subject to change)**

The following articles are affected by this change:

Cryo.s with Datamatrix Code and Linear Barcode

2 ml versions with internal thread

Item No.	Description
122263-2DG	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, NATURAL SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122263-2D3	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, NATURAL SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, CUSTOMISED, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122277-2DG	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, GREEN SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122277-2D3	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, GREEN SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, CUSTOMISED, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122278-2DG	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, YELLOW SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122278-2D3	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, YELLOW SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, CUSTOMISED, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122279-2DG	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, BLUE SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122279-2D3	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, BLUE SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, CUSTOMISED, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122280-2DG	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, RED SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128
122280-2D3	CRYO.S, 2 ML, PP, ROUND BOTTOM, INTERNAL THREAD, RED SCREW CAP, STARFOOT, NATURAL, STERILE, 100 PCS./BAG, CUSTOMISED, DATAMATRIX ECC 200, 14X14 DOTS, BARCODE 128

Please note that the guaranteed product properties, quality certificates and product dimensions will remain unchanged.

If you would like to discuss this product modification further, please contact your local Account Manager or email info.uk@gbo.com

Your sincerely
Greiner Bio-One