

Frickenhausen, 20th April 2022

INFOMAILING

DESIGN MODIFICATION MULTIPURPOSE BEAKER

Dear Sir or Madam,

We are continuously striving to optimise our product portfolio and to improve existing products. In this context, we would like to inform you about a design change of the multipurpose beakers.

To handle the changeover in the best possible way, the changed multipurpose beakers will have changed new item numbers.

After the current stock has been sold and starting with **end of May 2022** the multipurpose beakers have the following properties:

- Working volume: 120 ml
- Shelf-life: 5 years
- Sterile versions: sterility level SAL 10⁻⁶
- Additional multipurpose beaker with tamper-proof label (Item No.724413)

The previous integrity seal ring will be discontinued. Instead, there will be an additional multipurpose beaker with a tamper-proof label.



Current Design



Future Design

In particular, the following articles are affected by this change:

Multipurpose Beakers

PREVIOUS ITEM NUMBERS		NEW ITEM NUMBERS	
Item No.	Product Description	Item No.	Product Description
724320	URINE BEAKER 100 ML, YELLOW LID, WITH RING SCREW CLOSURE, BULK-PACKED, STERILE	724410	MULTIPURPOSE BEAKER, PP, 120 ML, GRADUATED, YELLOW LID, BULK-PACKED, STERILE, 300 PCS/CARTON
724402	MULTIPURPOSE BEAKER, PP, 100 ML, 61 X 71 MM, WITH SCREW CAP, GRADUATED, NATURAL, NON-STERILE, 300 PCS/BOX	724411	MULTIPURPOSE BEAKER, PP, 120 ML, GRADUATED, WHITE LID, BULK-PACKED, NON-STERILE, 300 PCS/CARTON
724461	MULTIPURPOSE BEAKER, PP, 100 ML, 61 X 71 MM, WITH SCREW CAP, GRADUATED, NATURAL, STERILE, 300 PCS/BOX	724412	MULTIPURPOSE BEAKER, PP, 120 ML, GRADUATED, WHITE LID, BULK-PACKED, STERILE, 300 PCS/CARTON
724473	MULTIPURPOSE BEAKER, PP, 100 ML, 61 X 71 MM, WITH SCREW CAP, GRADUATED, RSC-RED, STERILE, 300 PCS/BOX	724414	MULTIPURPOSE BEAKER, PP, 120 ML, GRADUATED, RED LID, BULK-PACKED, STERILE, 300 PCS/CARTON
-	-	724413	MULTIPURPOSE BEAKER, PP, 120 ML, GRADUATED, WHITE LID, BULK-PACKED, WITH TAMPER-PROOF LABEL , STERILE, 300 PCS/CARTON
724170	SCREW CAP, 61 MM, NATURAL, FOR 7244XX, 720 PCS/BOX		Immediate discontinuation: no substitute, please convert to assembled new product
724173	SCREW CAP, 61 MM, RED, FOR 724401, 720 PCS./BOX		
724401	MULTIPURPOSE BEAKER BOTTOM, PP,100 ML, 61 X 71 MM, WITHOUT SCREW CAP, GRADUATED, NATURAL, NON-STERILE, 720 PCS/BOX		

Should you have any questions, please do not hesitate to contact your export manager, local subsidiary or GBO distributor.

Kind regards,



p.p. Michael Wirth
Director Marketing & Product Management