

Stonehouse, March 2022


## INFOMAILING

# PRODUCTION MODIFICATION SCREW CAPS FOR CELLMASTER POLYSTYRENE ROLLER BOTTLES


Dear Customer,

We would like to inform you about an upcoming change regarding the standard and filter screw caps for roller bottles. On a lot-controlled basis starting from April 2022 there will be a small change to the screw cap width as shown below (detail A).


This is accompanied by a slight change in one drawing dimension (detail A). Please find the corresponding product drawings below.




Current standard screw cap



New standard screw cap



Current filter screw cap



New filter screw cap

A list of affected roller bottle item numbers can be found on the following page.

Please note that the **guaranteed product properties** for screw caps and roller bottles including the quality certificates and data sheets **remain unchanged**. In particular, the **HDPE raw material** of the caps **will not be changed**.

If you have any queries regarding this change, please contact your local Greiner Bio-One Account Manager or email [info.uk@gbo.com](mailto:info.uk@gbo.com)

Yours sincerely  
Greiner Bio-One.

## THE FOLLOWING CATALOGUE ITEMS ARE AFFECTED BY THIS CHANGE:

### Roller bottles

### Polystyrene

Ø: 122 mm, Graduation: yes, Raw material: PS, Surface treatment: TC, Sterile: +

Item no.	Flask design	Height	Growth area	Cap colour	Surface	Total volume	Cap design	Sterile	Qty. inner / outer
680660	1 X	271 mm	850 cm <sup>2</sup>	● blue	smooth	2,520 ml	screw cap	+	2 / 24
680665	1 X	271 mm	850 cm <sup>2</sup>	● blue	smooth	2,520 ml	screw cap	+	24 / 24
680648	1 X	271 mm	850 cm <sup>2</sup>	● blue	smooth	2,520 ml	screw cap	+	24 / 48
680658	1 X	271 mm	850 cm <sup>2</sup>	● blue	smooth	2,520 ml	filter screw cap	+	2 / 24
680668	1 X	271 mm	850 cm <sup>2</sup>	● blue	smooth	2,520 ml	filter screw cap	+	24 / 24
680645	1 X	271 mm	850 cm <sup>2</sup>	● blue	smooth	2,520 ml	filter screw cap	+	24 / 48
681670	2.5 X	271 mm	2,125 cm <sup>2</sup>	● blue	ribbed	2,300 ml	screw cap	+	2 / 24
681675	2.5 X	271 mm	2,125 cm <sup>2</sup>	● blue	ribbed	2,300 ml	screw cap	+	24 / 24
681672	2.5 X	271 mm	2,125 cm <sup>2</sup>	● blue	ribbed	2,300 ml	filter screw cap	+	2 / 24
682660	1 XL	500 mm	1,700 cm <sup>2</sup>	● blue	smooth	4,970 ml	screw cap	+	1 / 12
682612	1 XL	500 mm	1,700 cm <sup>2</sup>	● blue	smooth	4,970 ml	screw cap	+	12 / 12
682624	1 XL	500 mm	1,700 cm <sup>2</sup>	● blue	smooth	4,970 ml	screw cap	+	12 / 24
682625	1 XL	500 mm	1,700 cm <sup>2</sup>	● blue	smooth	4,970 ml	filter screw cap	+	12 / 24
682615	1 XL	500 mm	1,700 cm <sup>2</sup>	● blue	smooth	4,970 ml	filter screw cap	+	12 / 12
682670	5 XL	500 mm	4,250 cm <sup>2</sup>	● blue	ribbed	4,640 ml	screw cap	+	1 / 12
682672	5 XL	500 mm	4,250 cm <sup>2</sup>	● blue	ribbed	4,640 ml	screw cap	+	12 / 24
682678	5 XL	500 mm	4,250 cm <sup>2</sup>	● blue	ribbed	4,640 ml	filter screw cap	+	1 / 12

/ Sterile, individually packed screw caps: standard (item no. 383361) and filter (item no. 383382).