

Stonehouse, March 2022

# PRODUCT MODIFICATION

## UDI-CODE CHANGE


Dear Sir or Madam,


We would like to inform you that we are changing the **Unique Device Identification (UDI)-Code** for some of our products by adding the **date of manufacture**. In addition, there is a **new symbol** on the labels placed next to the numeric written date of manufacture. As a result, both the labels and the packaging layout have been adapted. The item numbers as well as the corresponding LOT-numbers indicating when the change will take effect can be found in the following table:



Product group	Item No.	LOT No.
VACUETTE® Multiple Use Drawing Needle	450076	22A04D
Safety Blood Collection Sets	450081	22A04
MiniCollect® Safety Lancet	450427	N222011
MiniCollect® Lancelino Safety Lancet	450520	B50G823G7


### OLD INNER BOX LABEL



**SAFETY Blood Collection Set + Holder**  
25G x 3/4", tubing length 12" (30 cm)

REF 450080  0123

LOT A0000000  36°C 97°F

 yyyy-mm-dd  4°C 40°F


GS1 UDI  (01)000000000000000  
(17)000000  
(10)A000000K


 STERILE  EO


### NEW INNER BOX LABEL

**SAFETY Blood Collection Set + Holder**  
25G x 3/4", tubing length 12" (30 cm)

REF 450080 24 pcs.

GS1 UDI  (01)000000000000000  
(11)000000  
(17)000000  
(10)A000000K

LOT 000000  0000-00-00

 0000-00-00

OLD CARTON LABEL

**VACUETTE®**

**SAFETY Blood Collection Set + Holder**

25G x 3/4", tubing length 12" (30 cm)  
single-packed, sterile, not made with natural rubber latex  
10 boxes of 24 pcs., 240 pcs. in total

REF **450080** CODE 128 (ALL CHARS)

LOT **A0000000** CODE 128 (ALL CHARS)

yyyy-mm-dd CODE 128 (ALL CHARS)

CE 0123 STERILE EO

Production location:  
Nipro (Thailand) Corporation Ltd., 10/2 Moo 8 Bangnomko, Sena, Phra Nakhon Si Ayutthaya 13110, Thailand  
Made in Thailand distributed by Greiner Bio-One GmbH, Austria

4°C 40°F 36°C 97°F

GS1 UDI

GS1 UDI Code (ALL CHARS)

(01)00000000000000(17)000000(10)A000000K

Greiner Bio-One GmbH  
Bad Haller Str. 32  
4550 Kremsmünster  
Austria

CAUTION: U. S. Law restricts this device to sale by or on the order of a physician.  
Medical device for venous blood collection and in vitro diagnostic use.

HIBC DATA-MATRIX

NEW CARTON LABEL

**VACUETTE®**

**SAFETY Blood Collection Set + Holder**

25G x 3/4", tubing length 12" (30 cm)  
single-packed, sterile, not made with natural rubber latex  
10 boxes of 24 pcs., 240 pcs. in total

REF **450080**

LOT **A0000000**

yyyy-mm-dd

yyyy-mm-dd

CE 0123 STERILE EO

Production location:  
Nipro (Thailand) Corporation Ltd., 10/2 Moo 8 Bangnomko, Sena, Phra Nakhon Si Ayutthaya 13110, Thailand  
Made in Thailand distributed by Greiner Bio-One GmbH, Austria

4°C 40°F 36°C 97°F

GS1 UDI

DATA-MATRIX (01)00000000000000 (11)YYMMDD (17)YYMMDD (10)N0000000

Greiner Bio-One GmbH  
Bad Haller Str. 32  
4550 Kremsmünster  
Austria

C/No.1

CAUTION: U. S. Law restricts this device to sale by or on the order of a physician.  
Medical device for venous blood collection and in vitro diagnostic use.

HIBC DATA-MATRIX

If you have any questions regarding this change, please contact your local Greiner Bio-One Account Manager or email [info.uk@gbo.com](mailto:info.uk@gbo.com)

Yours sincerely,  
Greiner Bio-One