



greiner bio-one

# CELLview™ Slide

Product of the Year 2017



## An exceptional product ...

- ... offering optimal test conditions
- ... with maximum optical transparency
- ... facilitating even ambitious experiments
- ... ensuring reliable research results



## ... that convinces!

- Innovative design
- Simple & efficient usage
- Perfectly suited for a variety of microscopic analyses
- Black compartmentalisation inhibits bleaching
- Reduced meniscus effect for optimal cell analyses

CELLview™ Slide is the ideal tool for cell culture, immunocytochemical analysis and high-resolution microscopy.

## At a glance

Black compartmentalisation



inhibits crosstalk and bleach out in neighbouring wells

96 well plate design



facilitates usage of multi-channel pipettes

Conical round wells



reduce meniscus effect

Embedded cover glass



guarantees an even focal plane

Surface treatment



plastic-like cell culture surfaces on glass for optimal cultivation of any type of cells

Alphanumeric well coding



ensures unmistakable sample identification

## Key Features

- 1** Detachable black compartmentalisation
- 2** 10 round wells equal to the design of a 96 well microplate
- 3** 0,17 mm high transparent borosilicate cover glass embedded in the slide
- 4** Surface treatment for optimal cell growth
- 5** Well coding on slide and compartmentalisation
- 6** Positioning notch for automated microscopy



For further information and/or sample ordering please visit our website [www.gbo.com](http://www.gbo.com) or contact us.

### Greiner Bio-One:

#### **Austria**

[office@at.gbo.com](mailto:office@at.gbo.com)

#### **Belgium**

[info@be.gbo.com](mailto:info@be.gbo.com)

#### **Brazil**

[office@br.gbo.com](mailto:office@br.gbo.com)

#### **China**

[office@cn.gbo.com](mailto:office@cn.gbo.com)

#### **France**

[accueil@gbo.com](mailto:accueil@gbo.com)

#### **Germany**

[info@de.gbo.com](mailto:info@de.gbo.com)

#### **Hungary**

[office@hu.gbo.com](mailto:office@hu.gbo.com)

#### **Italy**

[office@it.gbo.com](mailto:office@it.gbo.com)

#### **Japan**

[info@jp.gbo.com](mailto:info@jp.gbo.com)

#### **Netherlands**

[info@nl.gbo.com](mailto:info@nl.gbo.com)

#### **Portugal**

[info@vacuette.pt](mailto:info@vacuette.pt)

#### **Spain**

[info@vacuette.es](mailto:info@vacuette.es)

#### **UK**

[info@uk.gbo.com](mailto:info@uk.gbo.com)

#### **USA**

[office@us.gbo.com](mailto:office@us.gbo.com)



Contact information of Greiner Bio-One subsidiaries and distributors can be found on our website:  
[www.gbo.com/subsidiaries-distributors](http://www.gbo.com/subsidiaries-distributors)