



MiniCollect® Safety Lancets

Gentle & Caring for Patients Both Young and Old



- Safe and hygienic capillary blood collection
- More precision for simple sample collection
- Unique identification of different types with colour code
- Practical handling for increased collection comfort
- Can be applied for glucose testing

More Safety



- Absolute safety due to the automatic retraction of the blade / needle into the plastic housing after puncture
- Unintentional activation is practically impossible. The blade / needle is only activated when the trigger button is pushed.
- Designed for single-use only

Precision and Simple Handling



- Each version ensures a consistent puncture depth for straightforward sample collection
- The standardised blade / needle positioning means reliably consistent puncture alignment
- Ergonomic design for more handling comfort
- Self-explanatory application saves time
- Small, compact product design saves space in disposal container

Clear Identification

Available in five different puncture depths, which are identifiable with different colours.

Application



1 Remove the protective cover

2 Position the lancet

3 Puncture the skin by pressing button

4 Discard first drop of blood

Item Numbers

Item n°.	Description	Blood Volume	Inner pack	Outer pack
MiniCollect® Safety Lancets with blade				
450427	Safety Lancet, pink, penetration depth 1 mm	Low blood volume	200 pcs.	2000 pcs.
450428	Safety Lancet, green, penetration depth 1,5 mm	Medium blood volume	200 pcs.	2000 pcs.
450429	Safety Lancet, blue, penetration depth 2 mm	High blood volume	200 pcs.	2000 pcs.
MiniCollect® Safety Lancets with needle				
450438	Safety Lancet, lavender, 28 G, penetration depth 1,25 mm	Single blood drop	200 pcs.	2000 pcs.
450439	Safety Lancet, orange, 23 G, penetration depth 2,25 mm	Single blood drop	200 pcs.	2000 pcs.