



# VeinViewer® Flex

Illuminating Healthcare



- Universal - useful for all patients
- Adaptable to all situations
- Reliability in daily routine
- PDP: Pre-, During and Post-access advantages

## Features built for mobility



More and more patients with challenging venipuncture deserve a solution. VeinViewer® Flex is that solution to help meet those needs.

VeinViewer® Flex packs the latest advances in superior vascular imaging into a handheld, mobile package. VeinViewer® can help deliver faster venipuncture procedures on your difficult patients through greatly improved first stick success. Flex gives you revolutionary HD views in realtime, so you have the best view of vasculature procedures pre-, during and post-access. MaxDepth allows a view of peripheral veins up to 10mm deep and blood patterns up to 15mm beneath the surface of the skin.



### ✓ Fast Swap. Flexible Power.

Flex runs on Fast Swap Battery and AC power so there are no concerns about a down device.

### ✓ Hospital Grade Durability.

Durability withstands even the busiest shifts in the most demanding conditions.

### ✓ S-Mount (Stability Mount).

When you need to be more hands-on, the attachable S-mount gives you a hands-free approach.

### ✓ Flex Imaging Unit.

Built to be controlled and positioned easily with one hand, Flex is ergonomically designed for simplified procedures. The standard Universal Mode provides HD and Df<sup>2</sup> projection technology, a clear and accurate view of superficial vein patterns in real time, even in remote locations.

## ASSESS™ Imaging Suite

See veins like never before with the ASSESS™ Imaging Suite of software only available on VeinViewer® products. The choice of various modes ensures optimal visualization of veins for more control and better assessment of each individual patient.

### ☞ Universal

Useful for all patients

### ☞ Fine Detail

Critical for seeing small veins in pediatric or sclerotherapy patients

### ☞ Image Capture

HIPAA compliant, takes up to 200 static .PNG images of the blood pattern for export to EMR documentation

### ☞ Inverse

Reverses the colour of the veins within the image field for further image customization

### ☞ Resize

Choose 1 of 3 window sizes, a great feature for paediatric patients

### ☞ MaxBright

Brightest vein illumination image, needed for varying lighting

Item N°	Description	Packaging
890001	VeinViewer® Flex, Contains 2 pcs. Lithium Ion batteries 2 S-Mounts (Super Clamp, C-Clamp), 2 Fast Swap Batteries, Battery Charger, AC Power Adapter, International Power Adapters, Nylon Storage Case, USB Cable, Product Training DVD	all in one package
890010	Spare Rechargeable Battery for VeinViewer® Flex and Vision 2	-