



MiniCollect® Capillary Blood Collection Tubes

For your most precious samples



- Improved Cap Design
- Easier opening thanks to better grip
- Maximum tube opening for minimum adhesion
- Maximum insight due to transparent labels and colour-coded caps



MiniCollect® Capillary Blood Collection Tubes

Improved Cap Design

MiniCollect® offers a gentle way to collect small blood samples for a wide range of analyses. The system is perfectly suited for young children, geriatric patients as well as patients with fragile veins. **MiniCollect®** is also recommended in situations where venous blood collection proves particularly difficult, such as for patients with extensive burns.

Easier opening thanks to better grip

The **MiniCollect®** Capillary Blood Collection System is made for your most precious samples. The **MiniCollect®** Cap has been made easier to open.

The number of ridges has been increased and they have been made deeper to enable a better grip of the product. In addition, the cap rim overlaps the top of the tube to make it easier to grip and push off with the thumb.



More information at
shop.gbo.com