




PRODUCT MODIFICATION

New MiniCollect® Cap

Dear Customer,

As part of our ongoing optimisation of products we have introduced a new and improved **MiniCollect®** cap.

What's New?

-  An increased number of ridges
-  Deeper ridges to improve the **grip** of the cap
-  The cap rim overlaps the top of the tube allowing for an **easy push-off**

Please see the tables overleaf for a list of affected products:



Overview of affected MiniCollect® Tubes

Item No.	Description	Cap Colour	Fill Volume	NPC Code
450530	K3E K3EDTA	lavender	0.25 - 0.5 ml	KFK834
450531	K3E K3EDTA	lavender	1 ml	KFK835
450532	K2E K2EDTA	lavender	0.25 - 0.5 ml	KFK836
450542	Z No Additive	white	-	KFK798
450533	CAT Serum Separator	gold	0.5 - 0.8 ml	KFK793
450534	CAT Serum	red	0.5 - 1 ml	KFK828
450535	LH Lithium Heparin Separator	light green	0.8ml	KFK829
450536	LH Lithium Heparin	green	0.5 ml	-
450537	LH Lithium Heparin	green	1 ml	KFK827
450538	9NC Coagulation sodium citrate 3.2%	light blue	0.5 ml	KFK832
450539	9NC Coagulation sodium citrate 3.2%	light blue	1 ml	KFK833
450540	FX Sodium Fluoride / Potassium Oxalate	grey	0.25 ml	KFK830
450541	FX Sodium Fluoride / Potassium Oxalate	grey	0.5 ml	KFK831

Overview of affected MiniCollect® Complete Tubes

Item No.	Description	Cap Colour	Fill Volume	NPC Code
450555	Z No Additive	white	-	KFK811
450545	K3E K3EDTA	lavender	0.25 - 0.5 ml	KFK816
450546	K3E K3EDTA	lavender	1 ml	KFK817
450547	K2E K2EDTA	lavender	0.25 - 0.5 ml	KFK818
450548	CAT Serum Separator	gold	0.5 - 0.8 ml	KFK813
450549	CAT Serum	red	0.5 - 1 ml	KFK812
450550	LH Lithium Heparin Separator	light green	0.8 ml	KFK815
450551	LH Lithium Heparin	green	1 ml	KFK814
450552	FX Sodium Fluoride /Potassium Oxalate	grey	0.25 ml	KFK819

The new cap will be available for all MiniCollect® tubes mentioned in the tables above **starting with LOT A2002----**.

If you have any questions, please feel free to contact your Area Sales Manager.

Best regards
Greiner Bio-One Ltd.



Greiner Bio-One Ltd., Brunel Way,
Stroudwater Business Park, GL10 3SX
PHONE: 01453 825255
E-MAIL: info.uk@gbo.com