



greiner bio-one

# CONFIRMATION

---

**Product description: Cryo.s™, Cryo.s™ with Datamatrix, Cryo.s™ Biobanking Tubes**

**Item no: 12x2xx, 12x2xx-xxx, 12x2xx-TRI, 97x5xx, 97x5xx-xxx, 13x2xx, 13x2xx-xxx**


Our raw material supplier confirm:

„We inform you that traces of Di-isobutyl phthalate (DIBP), Dibutyl phthalate (DBP, as an impurity of DIBP) and/or Diethyl phthalate (DEP) may be present in some of our products, as catalyst modifier residues. The maximum residue would not exceed 1 ppm, an extremely low level compared the maximum level (1000ppm) authorized for substances mentioned in the Candidate List of Substances of Very High Concern (Regulation 1907/2006/EC – REACH). Please note that DIBP and DEP are neither included in nor targeted by the European Decision 2004/781/EC, concerning measures relating to toys prohibiting certain phthalates.”

This confirmation is valid until 26.04.2017.

26.04.2016

\_\_\_\_\_ date



\_\_\_\_\_  
Quality Assurance  
Greiner Bio-One GmbH  
Maybachstraße 2  
D-72636 Frickenhausen, Germany