



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	142C.USE
GENERAL DESCRIPTION	TRANS AMIES LQ A/ALL 500PKG
DETAILED DESCRIPTION	1 box containing 10 Vi-Pak, each containing 50 film-film pouches with a tube with 1.2 ml of liquid Amies Medium with white cap without receptacle and 1 aluminum swab with rayon tip and orange cap

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
AMIES LIQUID	Sodium chloride Potassium chloride Calcium Chloride Magnesium Chloride Monopotassium phosphate Disodium Phosphate sodium thioglycollate Distilled water

TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 154.50 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.2 mm EXTERNAL DIAMETER: 12.0 mm \pm 0.5 mm
SPONGE	SPONGE	DIMENSIONAL SPECIFICATIONS
		LENGTH: 37.50 mm + 2.0 mm -1.50 mm WIDTH: 10.0 mm \pm 0.5 mm THICKNESS: 5.0 mm \pm 0.5 mm

This bulletin (Rev.02-31/03/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

CAP	WHITE CAP W/O RECEPTACLE	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.0 mm \pm 0.10 mm CAP LENGTH: 32.80 mm \pm 0.20 mm INSERTION DIAMETER: 9.00 mm +0.05 mm
APPLICATOR	ALUMINIUM APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH (*purchased product, no tolerances): 149.0 mm STICK DIAMETER HANDLE PART: 0.9 mm \pm 0.1 mm TIP DIAMETER: 3.0 mm \pm 0.5 mm FIBRE TIP LENGTH: 15.0 mm \pm 2.0 mm TOTAL LENGTH: 150.0 mm \pm 1.0 mm
CAP	ORANGE CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.20 mm \pm 0.10 mm CAP LENGTH: 32.50 mm \pm 0.20 mm INSERTION DIAMETER: 9.30 mm + 0.05 mm

PACKAGING SPECIFICATION

OUTER BOX QTY	500
PACKAGING	FILM-FILM POUCH
VI-PAK QUANTITY	50

This bulletin (Rev.02-31/03/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET****IMPORTANT INFO**

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	15 MONTHS
MEDIUM pH	7.2 ± 0.3