


# Ordering Sheet Barcode Service for Cryo.s™

Orderer \_\_\_\_\_  
 Institution \_\_\_\_\_  
 Street \_\_\_\_\_  
 Postal code/City \_\_\_\_\_

Country \_\_\_\_\_  
 Phone \_\_\_\_\_  
 Fax \_\_\_\_\_  
 E-Mail \_\_\_\_\_

 For detailed information about "barcode" please ask for our brochure **Barcode-ABC** under (+49) 7022/948-0 or send an E-Mail to info@de.gbo.com.

Please tick off which kind of Cryo.s™ and screw cap colour you need:




Internal thread,  
1 ml volume,  
self standing




Internal thread,  
2 ml volume




Internal thread,  
2 ml volume,  
self standing




Internal thread,  
5 ml volume




External thread,  
2 ml volume,  
self standing




External thread,  
3.5 ml volume,  
self standing

Colour screw cap:

transparent

yellow

red

green

blue

- Selection of different barcode systems: Code 128; Code 39; 2/5 interleave; EAN; UPC and other
- Recommendation: Code 128C (only numerical code, represents two digits at a time)  
Advantage: compact data density and high reliability in reading
- New: Data Matrix (2D barcode) very high data density  
Advantage: possibility for clear identification (HIBC) of vial content
- All barcodes with check digit for internal verification and higher reliability in reading
- Printed barcode suitable for preservation in liquid nitrogen (-196 °C), resistant to ethanol (70 %) and DMSO (50 %)
- Minimum order: one original carton (300 pieces for large Cryo.s™ and 500 pieces for small Cryo.s™)

Barcode type: \_\_\_\_\_

Sequence start: \_\_\_\_\_

Sequence end: \_\_\_\_\_

Total number of digits: \_\_\_\_\_

Total amount of Cryo.s™ with barcodes: \_\_\_\_\_

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Released by (signature): \_\_\_\_\_

Date: \_\_\_\_\_

Greiner Bio-One GmbH  
 Maybachstraße 2  
 D-72636 Frickenhausen  
 Phone: (+49) 70 22/9 48-0  
 Fax: (+49) 70 22/9 48-514  
 E-Mail: info@de.gbo.com  
 www.gbo.com/bioscience

