

Cryo.s™ with Datamatrix

Selection of compatible scanners and others equipment from 3rd party suppliers

Task	Compatible Devices				
Full Rack Scanning, 48way rack	Ziath DataPaq™ Deep Focus Single Rack Scanner Version 2 Website	Ziath DataPaq™ Multirack Scanner Website	BioSero BioRead Deep Focus Single Rack Scanner Website	BioMicroLab SampleScan 96, Mini, HS and E-Series Website	Hamilton LabElite I.D. Reader Website Compatibility Statement
Full Rack Scanning, 81way rack	BioMicroLab SampleScan 96 Website	Ziath DataPaq™ Multirack Scanner Website			
Automated Decapping, 48way rack and internal thread tubes	Hamilton LabElite DeCapper Website Compatibility Statement	Hamilton LabElite I.D. Capper Website Compatibility Statement	Hamilton LabElite Integrated I.D. Capper Website Compatibility Statement		
Semi-automated and automated storage of 48way rack	Askion ASKION C-line® hermetic storage HS200 Website Compatibility Statement Video	Hamilton SAM, ASM, Verso, BIOS Website Compatibility Statement			
Liquid handling tubes handling and sample processing for 48way rack	Hamilton Microlab STAR and NIMBUS, VATNAGE Liquid Handling Workstations, Rack Runner Website Compatibility Statement				

Revision: June 2015