

# 96 Well ELISA Strip Plates


## For Diagnostics and Immunological Research

Greiner Bio-One is offering a large variety of 96 well strip plates for diagnostic and immunological research applications.

Microplates in strip format offer the advantage of greater flexibility in diagnostics. Individual strips can be removed from the support frame so that the number of tests to be performed can be adjusted to the number of samples, and is not predetermined by the microplate format used. In addition, the individual strips can be subjected to a wide variety of different test conditions. The quality of our immunological products is constantly controlled in our quality assurance laboratory by means of ELISA. Clear MICROLON® products are tested in a colorimetric immunoassay, black FLUOTRAC™ plates in a fluorescence immunoassay, white LUMITRAC™ plates in a luminescence immunoassay.

### Key Facts

- Strips mounted in frame
- U8 clear strips
- F8 clear / white / black strips
- F16 / U16 clear strips
- Black and white F16 / U16 strip plates on request
- L: 127.76 mm / W: 85.48 mm



Please turn page for ordering information!

Ordering Information			
Cat. No.	Product Description	Quantity per Bag	Quantity per Case
754 070	96 well strip plate, 6 x U16 strips, U-bottom, MICROLON® 200, med. binding, clear	10	100
754 061	96 well strip plate, 6 x U16 strips, U-bottom, MICROLON® 600, high binding, clear	10	100
756 070	96 well strip plate, 6 x F16 strips, F-bottom, MICROLON® 200, med. binding, clear	10	100
756 071	96 well strip plate, 6 x F16 strips, F-bottom, MICROLON® 600, high binding, clear	10	100
762 070	96 well strip plate, 12 x F8 strips, F-bottom, MICROLON® 200, med. binding, clear	10	100
762 071	96 well strip plate, 12 x F8 strips, F-bottom, MICROLON® 600, high binding, clear	10	100
762 075	96 well strip plate, 12 x F8 strips, F-bottom, LUMITRAC™ 200, med. binding, white	10	100
762 074	96 well strip plate, 12 x F8 strips, F-bottom, LUMITRAC™ 600, high binding, white	10	100
762 076	96 well strip plate, 12 x F8 strips, F-bottom, FLUOTRAC™ 200, med. binding, black	10	100
762 077	96 well strip plate, 12 x F8 strips, F-bottom, FLUOTRAC™ 600, high binding, black	10	100
767 070	96 well strip plate, 12 x U8 strips, U-bottom, MICROLON® 200, med. binding, clear	10	100
767 071	96 well strip plate, 12 x U8 strips, U-bottom, MICROLON® 600, high binding, clear	10	100