

Sapphire Pipette and Filter Tips

Maximum precision and recovery with liquid handling systems from Greiner Bio-One

The Sapphire product family comprises standard, filter and low retention tips manufactured from virgin materials. All tips are transparent, graduated, and allow precise and comfortable pipetting with maximal recovery.

Visual Identification System


Spot your tip instantly with Sapphire's new Visual Identification System: the boxes are transparent with coloured tip inserts based on the volume of the relevant tip. (Pipette tips have not been changed and still fit into traditional blue boxes.)

- 10 µl tips → Red insert
- 20 µl / 100 µl / 200 µl tips → Yellow insert
- 300 µl tips → Green insert
- 1250 µl tips → Blue insert


Key Facts

- Available in 10 µl / 20 µl / 100 µl / 200 µl / 300 µl / 1250 µl sizes
- Extended 10 µl tip for improved reach and better sample recovery
- Graduations for perfect visual control of liquid
- Thin-walled top for optimal seal
- Free of detectable DNase, RNase, human DNA and non-pyrogenic
- Low Retention Tips** for even higher recovery

Packaging Options	Rack	Refill	Bag
Pipette Tips	✓	✓	✓
Low Retention Pipette Tips	✓	✓	
Filter Tips	✓		✓
Low Retention Filter Tips	✓		✓



Sapphire Visual Identification System



Sapphire Pipette Tips

Ordering Information

Cat. No.		Product Description	Pieces/ Box	Pieces/ Case
Pipette Tips				
Rack	771 254	10 µl Pipette Tips, graduated, racked, PP, natural, non-sterile, 96 tips/rack	960	4800
	737 254	200 µl Pipette Tips, graduated racked, PP, natural, non-sterile, 96 tips/rack	960	4800
	738 254	300 µl Pipette Tips, graduated, racked, PP, natural, non-sterile, 96 tips/rack	960	4800
	750 254	1250 µl Pipette Tips, graduated, racked, PP, natural, non-sterile, 96 tips/rack	768	3840
Refill	771 257	10 µl Pipette Tips, graduated, refill, PP, natural, non-sterile, 96 tips/insert	960	4800
	737 257	200 µl Pipette Tips, graduated, refill, PP, natural, non-sterile, 96 tips/insert	960	4800
	738 257	300 µl Pipette Tips, graduated, refill, PP, natural, non-sterile, 96 tips/insert	960	4800
	750 257	1250 µl Pipette Tips, graduated, refill, PP, natural, non-sterile, 96 tips/insert	768	3840
Bag	771 250	10 µl Pipette Tips, graduated, bagged, PP, natural, non-sterile, 1000 tips/bag	1000	10000
	738 250	300 µl Pipette Tips, graduated, bagged, PP, natural, non-sterile, 1000 tips/bag	1000	10000
	750 250	1250 µl Pipette Tips, graduated, bagged, PP, natural, non-sterile, 768 tips/bag	768	7680
Low Retention Pipette Tips				
Rack	771 255	10 µl Low Retention Pipette Tips, grad., racked, PP, natural, non-sterile, 96 tips/rack	960	4800
	737 255	200 µl, Low Retention Pipette tips, grad. racked, PP, natural, non-sterile, 96 tips/rack	960	4800
	738 255	300 µl Low Retention Pipette tips, grad., racked, PP, natural, non-sterile, 96 tips/rack	960	4800
	750 255	1250 µl Low Retention Pipette tips, grad., racked, PP, natural, non-sterile, 96 tips/rack	768	3840
Refill	771 258	10 µl Low Retention Pipette Tips, grad., refill, PP, natural, non-sterile, 96 tips/insert	960	4800
	737 258	200 µl, Low Retention Pipette Tips, grad. refill, PP, natural, non-sterile, 96 tips/insert	960	4800
	738 258	300 µl Low Retention Pipette Tips, grad., refill, PP, natural, non-sterile, 96 tips/insert	960	4800
	750 258	1250 µl Low Retention Pipette Tips, grad., refill, PP, natural, non-sterile, 96 tips/insert	768	3840
Filter Tips				
Rack	771 261	10 µl Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	773 261	20 µl Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	737 261	100 µl, Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	738 261	300 µl Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	750 261	1250 µl Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	768	3840
Bag	771 260	10 µl Filter Tips, graduated, bagged, PP, natural, sterile, 1000 tips/bag	1000	10000
	738 260	300 µl Filter Tips, graduated, bagged, PP, natural, sterile, 1000 tips/bag	1000	10000
	750 260	1250 µl Filter Tips, graduated, bagged, PP, natural, sterile, 768 tips/bag	768	7680
Low Retention Filter Tips				
Rack	771 265	10 µl Low Retention Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	773 265	20 µl Low Retention Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	New 772 265	100 µl Low Retention Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	New 737 265	200 µl Low Retention Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	738 265	300 µl Low Retention Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	960	4800
	750 265	1250 µl Low Retention Filter Tips, graduated, racked, PP, natural, sterile, 96 tips/rack	768	3840
Bag	New 771 266	10 µl Low Retention Filter Tips, graduated, bagged, PP, natural, sterile, 1000 tips/bag	1000	10000
	New 737 266	200 µl Low Retention Filter Tips, graduated, bagged, PP, natural, sterile, 1000 tips/bag	1000	10000
	New 750 266	1250 µl Low Retention Filter Tips, graduated, bagged, PP, natural, sterile, 768 tips/bag	768	7680
Racks				
	970 310	Racks for 10 µl tips, unfilled, PP, natural with lid, red insert, for 96 tips	10	50
	970 320	Racks for 20 µl / 100 µl / 200 µl tips, unfilled, PP, natural with lid, yellow insert, for 96 tips	10	50
	970 330	Racks for 300 µl tips, unfilled, PP, natural with lid, green insert, for 96 tips	10	50
	970 350	Racks for 1250 µl tips, unfilled, PP, natural with lid, blue insert, for 96 tips	8	40