

# Centrifuge Tubes

## 15 ml / 50 ml polypropylene tubes

- ☞ Centrifuge tubes with conical bottom, graduation and writing area
- ☞ With screw cap
- ☞ Optimal mechanical, thermal and chemical stability
- ☞ Light Protection Tubes for light-sensitive materials and reactions

Maximum relative centrifugal forces\*:

	1882xx	2102xx	227261	227270	22728x
Swinging-bucket rotor	4000 x g	2800 x g	3200 x g	3200 x g	3200 x g
Fixed-angle rotor	15000 x g	9500 x g	17000 x g	17000 x g	9000 x g

\* The stated maximum RCF values are guidelines only, depending on a variety of factors such as rotor, temperature, density, pH and type of liquid. The suitability of tubes for a specific application using high centrifugation forces has to be tested.

### Key Facts

- ☞ Polypropylene
- ☞ Graduation and writing area
- ☞ Bulk or rack packaging available
- ☞ Brown Light Protection Tubes



Free of detectable  
DNase, RNase,  
human DNA  
non-pyrogenic

non-  
cytotoxic

## Ordering Information

Order No.	Product Description	Quantity per Bag	Quantity per Case
188261	15 ml tubes, sterile, graduation and writing area	50 (Rack)	500
188271	15 ml tubes, sterile, graduation and writing area	100	1000
188283	15 ml tubes, light protection tubes, sterile, graduation and writing area	50 (Rack)	500
188280	15 ml tubes, light protection tubes, sterile, graduation and writing area	100	1000
210261	50 ml tubes, with support skirt, sterile, graduation and writing area	25	450
210270	50 ml tubes, with support skirt, sterile, graduation and writing area	25 (Rack)	300
227261	50 ml tubes, sterile, graduation and writing area	20	500
227270	50 ml tubes, sterile, graduation and writing area	25 (Rack)	300
227285	50 ml tubes, sterile, graduation, writing area, white screw cap	20	500
227283	50 ml tubes, light protection tubes, sterile, graduation and writing area	25 (Rack)	300
227280	50 ml tubes, light protection tubes, sterile, graduation and writing area	20	500

For further information and/or sample ordering please visit our website [www.gbo.com](http://www.gbo.com) or contact us:

Germany (Main office): [info@de.gbo.com](mailto:info@de.gbo.com) | Austria: [office@at.gbo.com](mailto:office@at.gbo.com) | Belgium: [info.be@gbo.com](mailto:info.be@gbo.com) | Brazil: [office@br.gbo.com](mailto:office@br.gbo.com) | China: [info@cn.gbo.com](mailto:info@cn.gbo.com)  
 France: [accueil@gbo.com](mailto:accueil@gbo.com) | Hungary: [office@hu.gbo.com](mailto:office@hu.gbo.com) | Italy: [office@it.gbo.com](mailto:office@it.gbo.com) | Japan: [info.jp@gbo.com](mailto:info.jp@gbo.com) | Netherlands: [info.nl@gbo.com](mailto:info.nl@gbo.com)  
 Portugal: [info@vacuette.pt](mailto:info@vacuette.pt) | Spain: [info@es.gbo.com](mailto:info@es.gbo.com) | UK: [info.uk@gbo.com](mailto:info.uk@gbo.com) | USA: [office@us-gbo.com](mailto:office@us-gbo.com)

[www.gbo.com](http://www.gbo.com)