

Highly Transparent Centrifuge Tubes

12 and 15 ml Polystyrene Centrifuge Tubes

Besides the regular polypropylene tubes, Greiner Bio-One offers highly transparent 12 and 15 ml polystyrene tubes with screw caps.

- 12 ml tubes with conical bottom and supporting skirt
- 12 and 15 ml tubes with round bottom
- 15 ml centrifuge tubes with conical bottom, graduation and writing area

The tubes withstand the following centrifugation forces:

	163 1xx	164 1xx	186 1xx	188 1xx
ø (mm) x height (mm)	17 x 100	16.8 x 100	17 x 120	17 x 120
Swinging-bucket rotor	3,000 xg	5,000 xg	2,500 xg	2,500 xg
Fixed-angle rotor	5,000 xg	5,000 xg	2,800 xg	4,500 xg

Key Facts

- Highly transparent
- Graduation and writing area (only 15 ml tubes with conical bottom, Cat. No. 188 xxx)
- Suitable for centrifugation



Ordering Information

Cat. No.	Product Description	Quantity per Bag	Quantity per Case
163 177	12 ml centrifuge tube, round bottom, non-sterile, black cap	-	900
164 161	12 ml centrifuge tube, conical bottom, supporting skirt, sterile, blue cap	25	1000
164 180	12 ml centrifuge tube, conical bottom, supporting skirt, non-sterile, white cap	-	1300
186 161	15 ml tube, round bottom, sterile	50 (Box)	500
186 171	15 ml tube, round bottom, sterile	100	1000
188 161	15 ml centrifuge tube, conical bottom, graduation, writing area, sterile	50 (Box)	500
188 171	15 ml centrifuge tube, conical bottom, graduation, writing area, sterile	100	1000