

# Cryo Storage Boxes and Racks

## For Cryo.s™ Sample Tubes


The Cryo Storage Box is a helpful accessory for storing Cryo.s™ sample tubes at low and ultra-low temperatures. The box material polypropylene is very temperature and shock-resistant, thus allowing storage temperatures as low as -90 °C.

- ☞ Holds 81 Cryo.s™ sample tubes
- ☞ Translucent lid for optimum visibility of box content
- ☞ Venting holes for facilitated air exchange between interior and exterior
- ☞ Easy lid removal

The Cryo Storage Box fits all 1 ml and 2 ml Cryo.s™ with internal and external thread (Cat.No. 123 2XX, 121 2XX, 122 2XX und 126 2XX). It is compatible with standard metal racks.

### Key Facts

- ☞ Length: 126.5 mm
- ☞ Width: 126.5 mm
- ☞ Height\*: 51.0 mm
- ☞ For storage of 1 ml and 2 ml Cryo.s™ sample tubes
- ☞ 9 x 9 division to hold 81 tubes
- ☞ Temperature resistant polypropylene
- ☞ Temperature range: -90 °C ~ +121 °C



### Ordering Information

Cat. No.	Product Description	Quantity per Bag	Quantity per Case
802 202	Cryo Storage Box natural, lid natural	5	20
802 203	Cryo Storage Box red, lid natural	5	20
802 204	Cryo Storage Box blue, lid natural	5	20
802 206	Cryo Storage Box yellow, lid natural	5	20
802 225	Cryo Storage Box green, lid natural	5	20
802 501	Cryo One-Hand Operation Support Rack, polycarbonate, blue, with rubber base, for up to 40 Cryo.s™ with starfoot, Length: 200 mm, Width: 100 mm, Height: 22 mm	1	1

\* The indicated height refers to the Cryo Storage Box filled with 2 ml Cryo.s™ and covered with a lid.