

ThinCert™

Tissue Culture Inserts for Multiwell Plates

Greiner Bio.One offers tissue culture inserts for 6, 12 and 24 well multiwell plates with 0.4 µm, 1.0 µm, 3.0 µm and 8.0 µm pore sizes.

Advantages:

- improved cell adhesion due to physical surface treatment
- simplified pipetting due to self-lift geometry
- notches at the upper edge allow optimised gas exchange

Key Facts

- stable housing made of highly transparent polystyrene
- hanging geometry
- sealed PET capillary pore membrane
- single sterile blister packing



Free of detectable
DNase, RNase,
human DNA
non-pyrogenic

non-
cytotoxic

Ordering Information

Size	Cat. No.	Membrane Material	Pore Ø (µm)	Culture Surface (mm ²)	Pore Density (cm ⁻²)	Optical Membrane Properties	Plates/Box	Inserts/Box
6 Well	657 640	PET	0.4	452.4	1x10 ⁸	translucent	4	24
	657 641	PET	0.4	452.4	2x10 ⁶	transparent	4	24
	657 610	PET	1.0	452.4	2x10 ⁶	transparent	4	24
	657 630	PET	3.0	452.4	0.6x10 ⁶	transparent	4	24
	657 631	PET	8.0	452.4	2x10 ⁶	translucent	4	24
	657 638	PET	8.0	452.4	0.15x10 ⁶	translucent	4	24
12 Well	665 640	PET	0.4	113.1	1x10 ⁸	translucent	4	48
	665 641	PET	0.4	113.1	2x10 ⁶	transparent	4	48
	665 610	PET	1.0	113.1	2x10 ⁶	transparent	4	48
	665 630	PET	3.0	113.1	0.6x10 ⁶	transparent	4	48
	665 631	PET	8.0	113.1	2x10 ⁶	translucent	4	48
	665 638	PET	8.0	113.1	0.15x10 ⁶	translucent	4	48
24 Well	662 640	PET	0.4	33.6	1x10 ⁸	translucent	2	48
	662 641	PET	0.4	33.6	2x10 ⁶	transparent	2	48
	662 610	PET	1.0	33.6	2x10 ⁶	transparent	2	48
	662 630	PET	3.0	33.6	0.6x10 ⁶	transparent	2	48
	662 631	PET	8.0	33.6	2x10 ⁶	translucent	2	48
	662 638	PET	8.0	33.6	0.15x10 ⁶	translucent	2	48

For further information and/or sample ordering please visit our website www.gbo.com or contact us.

Germany (Main office): info@de.gbo.com | Austria: office@at.gbo.com | Belgium: info@be.gbo.com | Brazil: office@br.gbo.com | China: office@cn.gbo.com
 France: accueil@gbo.com | Hungary: office@hu.gbo.com | Italy: office@it.gbo.com | Japan: info@jp.gbo.com | Netherlands: info@nl.gbo.com
 Portugal: info@vacuette.pt | Spain: info@vacuette.es | UK: info@uk.gbo.com | USA: info@us.gbo.com

www.gbo.com/bioscience