

# 96 Well ELISA Microplates

## For Diagnostics and Immunological Research


Greiner Bio-One is offering a large variety of 96 well microplates for diagnostic and immunological research applications

- ☞ with U-Bottom / V-Bottom / F-Bottom standard (ST) / F-bottom chimney well / half area well geo-metry
- ☞ with medium and high binding surfaces
- ☞ made of clear / black / white polystyrene for transmission, fluorescence or luminescence measurements.

The quality of our immunological products is constantly controlled in our quality assurance laboratory by means of ELISA. Clear MICROLON® products are tested in a colorimetric immunoassay, black FLUOTRAC™ plates in a fluorescence immunoassay, white LUMITRAC™ plates in a luminescence immunoassay.

### Key Facts

- ☞ 96 well format
- ☞ U-/V-/ F-Bottom standard and chimney well
- ☞ Half Area well geometry
- ☞ Medium / high binding surfaces
- ☞ Clear / black / white
- ☞ L: 127.76 mm / W: 85.48 mm



Ordering information on page 2!

Ordering Information			
Cat. No.	Product Description	Quantity per Bag	Quantity per Case
650 001	Microplate, 96 well, PS, U-bottom, MICROLON® 200, med. binding, clear	10	40
650 061	Microplate, 96 well, PS, U-bottom, MICROLON® 600, high binding, clear	10	40
651 001	Microplate, 96 well, PS, V-bottom, MICROLON® 200, med. binding, clear	10	40
651 061	Microplate, 96 well, PS, V-bottom, MICROLON® 600, high binding, clear	10	40
655 001	Microplate, 96 well, PS, F-bottom/ST, MICROLON® 200, med. binding, clear	10	40
655 061	Microplate, 96 well, PS, F-bottom/ST, MICROLON® 600, high binding, clear	10	40
655 075	Microplate, 96 well, PS, F-bottom/chimney well, LUMITRAC™ 200, med. binding, white	10	40
655 074	Microplate, 96 well, PS, F-bottom/chimney well, LUMITRAC™ 600, high binding, white	10	40
655 076	Microplate, 96 well, PS, F-bottom/chimney well, FLUOTRAC™ 200, med. binding, black	10	40
655 077	Microplate, 96 well, PS, F-bottom/chimney well, FLUOTRAC™ 600, high binding, black	10	40
655 080	Microplate, 96 well, PS, F-bottom/chimney well, MICROLON® 200, med. binding, clear	10	40
655 081	Microplate, 96 well, PS, F-bottom/chimney well, MICROLON® 600, high binding, clear	10	40
675 001	Microplate, 96 well, PS, half area, MICROLON® 200 med. binding, clear	10	40
675 061	Microplate, 96 well, PS, half area, MICROLON® 600 high binding, clear	10	40