


Racks for Pipette Tips


Racks for Sapphire Tips (Bulk / Refill Units)



Cat. No.: 970 310



Cat. No.: 970 320



Cat. No.: 970 330



Cat. No.: 970 350

Racks for the EasyLoad® Refill System



Cat. No.: 941 300



Cat. No.: 941 310



Cat. No.: 941 320

Racks for Standard Pipette Tips (Bulk)



Cat. No.: 941 305



Cat. No.: 973 276



Cat. No.: 973 272



Cat. No.: 973 202



Cat. No.: 973 270



Cat. No.: 941 315



Cat. No.: 974 280



Cat. No.: 974 290



Cat. No.: 941 325

For further information and/or sample ordering please visit our website www.gbo.com or contact us.

Germany (Main office): info@de.gbo.com | Austria: office@at.gbo.com | Belgium: info@be.gbo.com | Brazil: office@br.gbo.com | China: office@cn.gbo.com
France: accueil@gbo.com | Hungary: office@hu.gbo.com | Italy: office@it.gbo.com | Japan: info@jp.gbo.com | Netherlands: info@nl.gbo.com
Portugal: info@vacuette.pt | Spain: info@vacuette.es | UK: info@uk.gbo.com | USA: info@us.gbo.com

Compatibility Pipette Tips - Empty Racks

Volume	Cat. No.	Sapphire Racks				EasyLoad® Racks				Standard Racks							
		970 310	970 320	970 330	970 350	941 300	941 310	941 320	941 305	941 315	941 325	973 202	973 270	973 272	973 276	974 280	974 290
Sapphire Pipette and Filter Tips																	
10 µl	771 250	•															
	771 257	•															
	771 258	•															
	771 260	•															
200 µl	737 257		•														
	737 258		•														
300 µl	738 250			•													
	738 257			•													
	738 258			•													
	738 260			•													
1250 µl	750 250				•												
	750 257				•												
	750 258				•												
	750 260				•												
EasyLoad® Pipette Tips																	
10 µl	741 015					•											
	741 020					•											
200 µl	741 000						•										
	741 061						•										
	741 065						•										
	741 070						•										
100-1000 µl	741 035							•									
	741 040							•									
	741 045							•									
	741 050							•									
Standard Pipette Tips																	
0,5-20 µl	765 290													•			
	771 290														•		
	771 291							•									
200 µl	685 290											•					
	685 295											•					
	739 290												•				
	739 291													•			
	739 296														•		
1000 µl	686 290																•
	686 295																•
	740 290																•
	740 291																•
	740 296																•
Conductive Pipette Tips																	
	686 296																•

For further information and/or sample ordering please visit our website www.gbo.com or contact us.

Germany (Main office): info@de.gbo.com | Austria: office@at.gbo.com | Belgium: info@be.gbo.com | Brazil: office@br.gbo.com | China: office@cn.gbo.com
 France: accueil@gbo.com | Hungary: office@hu.gbo.com | Italy: office@it.gbo.com | Japan: info@jp.gbo.com | Netherlands: info@nl.gbo.com
 Portugal: info@vacuette.pt | Spain: info@vacuette.es | UK: info@uk.gbo.com | USA: info@us.gbo.com