

CELLSTAR® Serological Pipettes

Serological Pipettes 1 to 50 ml

Greiner Bio-One offers in its CELLSTAR® product line a wide range of different serological pipettes. The volume capacity of the pipettes is increased by a negative graduation. The pipettes are provided with colour codes according to international standards. The use of high-grade polystyrene guarantees maximum clarity. Expiry date and LOT number are printed on each single packaging unit. All pipettes are supplied with a filter for protection against the suction of liquid into the pipetting device

- Design of the pipette tips guarantees a drop-free pipetting
- Maximum accuracy
- Clear, easily legible graduation
- A coloured vertical Schellbach stripe on 5, 10 and 25 ml pipettes makes it considerably easier to read off the volume

Key Facts

- Sterile
- High optical clarity
- 1 ml, 2 ml, 5 ml, 10 ml, 25 ml and 50 ml pipettes available
- Obtainable as bulk packaging, paper/plastic packaging and plastic/plastic packaging



Free of detectable
DNase, RNase,
human DNA
non-pyrogenic

non-
cytotoxic

Ordering Information

Cat. No.	Product Description	Quantity per Bag	Quantity per Case
604 107	1 ml pipette, 1/100 graduation, bulk packaging, sterile	25	1000
604 181	1 ml pipette, 1/100 graduation, paper/plastic packaging, sterile	1	1000
604 160	1 ml pipette, 1/100 graduation, plastic/plastic packaging, sterile	1	1000
710 107	2 ml pipette, 1/100 graduation, bulk packaging, sterile	25	1000
710 180	2 ml pipette, 1/100 graduation, paper/plastic packaging, sterile	1	1000
710 160	2 ml pipette, 1/100 graduation, plastic/plastic packaging, sterile	1	1000
606 107	5 ml pipette, 1/10 graduation, bulk packaging, sterile	25	500
606 180	5 ml pipette, 1/10 graduation, paper/plastic packaging, sterile	1	200
606 160	5 ml pipette, 1/10 graduation, plastic/plastic packaging, sterile	1	200
607 107	10 ml pipette, 1/10 graduation, bulk packaging, sterile	25	500
607 180	10 ml pipette, 1/10 graduation, paper/plastic packaging, sterile	1	200
607 160	10 ml pipette, 1/10 graduation, plastic/plastic packaging, sterile	1	200
760 107	25 ml pipette, 2/10 graduation, bulk packaging, sterile	25	200
760 180	25 ml pipette, 2/10 graduation, paper/plastic packaging, sterile	1	200
760 160	25 ml pipette, 2/10 graduation, plastic/plastic packaging, sterile	1	200
768 180	50 ml pipette, 1/2 graduation, paper/plastic packaging, sterile	1	100
768 160	50 ml pipette, 1/2 graduation, plastic/plastic packaging, sterile	1	100
	Shorties and Special Models		
710 183	2 ml aspiration pipette, no graduation, paper/plastic packaging, sterile, Special: no plug	1	1000
606 190	5 ml pipette, 1/10 graduation, plastic packaging, sterile, Special: shorty	1	200
607 190	10 ml pipette, 2/10 graduation, plastic packaging, sterile, Special: shorty	1	200