

# Disposal Bags

## Hygienic disposal of contaminated laboratory items

The disposal bags made of polypropylene (PP) foil are suitable for the sterilisation in steam autoclaves. The disposal bags made of polyamide (PA) are suitable for the hot-air sterilisation up to 160 °C. On request, polypropylene disposal bags are also available with imprint "Biohazard" if a sufficient number is ordered.

### Key Facts

- Available in polypropylene and polyamide
- Foil thickness: 0.05 mm
- Suitable for steam autoclaves (PP)
- Suitable for hot-air sterilisers (PA)
- Nominal capacities from 10 l – 130 l



### Ordering Information

Cat. No.	Product Description	Quantity per Case
643 201	Disposal bags, PP, 300 x 500 mm; nominal capacity 10 l	500
644 201	Disposal bags, PP, 400 x 780 mm; nominal capacity 30 l	500
646 201	Disposal bags, PP, 600 x 780 mm; nominal capacity 65 l	500
649 201	Disposal bags, PP, 700 x 1100 mm; nominal capacity 130 l	350
643 401	Disposal bags, PA, 300 x 500 mm; nominal capacity 10 l	500
644 401	Disposal bags, PA, 400 x 780 mm; nominal capacity 30 l	500
646 401	Disposal bags, PA, 600 x 780 mm; nominal capacity 65 l	300
649 401	Disposal bags, PA, 700 x 1100 mm; nominal capacity 130 l	200