


Ultra Low Profile Lid

Universal Ultra Low Profile Lid for LoBase Microplates

The ultra low profile lid without recesses in the edge of the lid are compatible with a wide range of different microplates. Manufactured of clear polystyrene, the ultra low profile lid is specifically recommended to seal 384 well Small Volume™ LoBase or 1536 well microplates.

With exceptional tolerances for suitability in automated systems, this lid guarantees the greatest possible protection against evaporation. Recesses in the edge of the lid improve readability of barcodes, and sealed microplates can easily be gripped from the side for transportation without problems.

Ultra low profile lids are available sterile or non-sterile, and customised packaging is available upon request.

<p>Key Facts</p> <ul style="list-style-type: none"> ➤ Ultra low profile lid without recesses ➤ Suitable for all microplates ➤ Robot-friendly ➤ Improved access for barcode reading ➤ Reduced evaporation ➤ Available sterile and non-sterile ➤ Length x Width x Height: 127.4 x 85.1 x 5.6 mm 	
---	--

Ordering Information			
Cat. No.	Product Description	Quantity per Bag	Quantity per Case
691 101	Universal Ultra Low Profile Lid, polystyrene, clear, non-sterile	5	100
691 161	Universal Ultra Low Profile Lid, polystyrene, clear, sterile	5	100