

Media Bottle

Tetragonal media bottle made of ultra-clear polyethylene terephthalate (PET)

Media bottle with screw cap for liquid storage. Thermal stability from -60 °C up to +60 °C.

Available in three sizes: 100, 500 and 1000 ml:

Key Facts	
<ul style="list-style-type: none"> ➤ 100 ml flask: 108.3 x 53.0 mm (height x width) ➤ 500 ml Flasche: 177.3 x 77.0 mm (height x width) ➤ 1000 ml Flasche: 218.8 x 92.0 mm (height x width) ➤ With graduation ➤ Sterile ➤ Triple-packed for GMP compliant workflow 	

Free of detectable
DNase, RNase,
human DNA
non-pyrogenic

non-
cytotoxic

Ordering Information

Cat. No.	Product Description	Quantity per Case
951 700	Media Bottle 100 ml	100
950 700	Media Bottle 500 ml	50
952 700	Media Bottle 1000 ml	24