

CX™ Extraction Kit

The CX™ Extraction Kit is designed to be used for rapid automated DNA preparation from cervical smears in combination with a platform of the Greiner Bio-One CX™ series (CX™ NIMBUS®, CX™ NIMBUS® plus, CX™ STARlet, and CX™ STARlet plus). It allows the preparation of high-quality DNA suitable for the direct downstream analysis by the PapilloCheck® (Greiner Bio-One; Cat. No. 465 060).

The CX™ Extraction Kit allows the automated fast and efficient purification of a multiple of 48 samples.

Key Facts

- Ready-to-use kit for 48 DNA preparations
- Preparation of high-quality DNA
- Validated in combination with PapilloCheck®
- Optimised for the DNA extraction from most common cervical sample collection systems



Application for PapilloCheck® (Cat. No. 465 060)

The CX™ Extraction Kit is validated for the isolation of viral and human DNA from human cervical smears collected with one of the following collection systems or media:

- PapilloCheck® Collection Kit, (Greiner Bio-One; Cat. No. 465 075)
- PreservCyt® (Hologic, Bedford, MA, USA)
- Surepath™ (BD, Franklin Lakes, NJ, USA)
- STM™ (Qiagen, Hilden, Germany)

Ordering Information

Cat. No.	Product Description	Content
517 070	CX™ Extraction Kit	48 preps

For further information and/or sample ordering please visit our website www.gbo.com or contact us.

Germany (Main office): info@de.gbo.com | Austria: office@at.gbo.com | Belgium: info@be.gbo.com | Brazil: office@br.gbo.com | China: office@cn.gbo.com
 France: accueil@gbo.com | Hungary: office@hu.gbo.com | Italy: office@it.gbo.com | Japan: info@jp.gbo.com | Netherlands: info@nl.gbo.com
 Portugal: info@vacuette.pt | Spain: info@vacuette.es | UK: info@uk.gbo.com | USA: info@us.gbo.com