



**"ACT"IONS
SPEAK LOUDER
THAN WORDS**



CELLSTAR®

CENTRIFUGE TUBES

**15 ml and 50 ml centrifuge tubes with the ACT® label -
Make an eco-conscious choice**

By choosing our ACT® certified centrifuge tubes, you are not only opting for excellent quality but also taking a significant step towards transparency and environmental responsibility with us.

The ACT® Label is independently audited and emphasizes **Accountability, Consistency and Transparency** to help you to make eco-conscious decisions that align with your sustainability goals.

STERILE	FREE OF detectable DNase	FREE OF detectable human DNA
FREE OF detectable RNase	non-cytotoxic	non-pyrogenic


KEY FACTS:

- / 15 ml and 50 ml tubes for universal applicability
- / Excellent thermal, mechanical and chemical stability
- / **Max. applicable relative centrifugation force of 15,000 x g (15 ml) and 17,000 x g (50 ml)**
(see product data sheet for details)
- / **Cartons made of ≥99 % recycled material**
- / **Made in Germany**
- / **15 ml:**
188261-N: 50 pcs. per rack
188271-N: 50 pcs. per bag
- / **50 ml:**
227261-N: 25 pcs per bag
227270-N: 25 pcs per rack
(available soon)



EXPLANATION OF THE ACT® LABEL

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERING PLEASE VISIT OUR WEBSITE OR CONTACT US.



Regional labels capture the differences in the shipping and end of life impacts for each region (US, EU, UK) in which the product is sold.

Manufacturing Impact Reduction
Reduced energy usage and waste generation

Renewable Energy Use
Production site has installed solar panels and participates in its utility provider's renewable energy program.

Shipping Impact
Made in Europe – optimized delivery routes

Packaging Content
Cartons made of ≥99 % recycled material

Lowest Environmental Impact Factor for 15 ml and 50 ml centrifuge tubes in Europe

PRODUCT OVERVIEW

Centrifuge Tubes

Polypropylene – 15 ml

Height: 120 mm, ø: 17 mm, Support skirt: no, Graduation: yes, Writing field: yes, Raw material: PP, Bottom Shape: conical, Working volume: ≤14 ml, Nominal volume: 15 ml, Cap design: screw cap, Sterile: +

Item No.	Product colour	Cap colour	Colour writing field	Graduation colour	Type of packaging	Sterile	Qty. inner / outer
188261-N	○ natural	● blue	○ white	● blue	rack	+	50 / 500
188271-N	○ natural	● blue	○ white	● blue	bag	+	50 / 1.000

Centrifuge Tubes

Polypropylene – 50 ml

Height: 115 mm, ø: 30 mm, Graduation: yes, Writing field: yes, Raw material: PP, Bottom Shape: conical, Working volume: ≤50 ml, Nominal volume: 50 ml, Cap design: screw cap, Sterile: +

Item No.	Product colour	Cap colour	Colour writing field	Graduation colour	Type of packaging	Sterile	Qty. inner / outer
227261-N	○ natural	● blue	○ white	● blue	bag	+	25 / 500
227270-N	○ natural	● blue	○ white	● blue	rack	+	25 / 250